



SWAB IDENTIFICATION CHART



1 Amies swab for bacterial culture



2 Aptima vaginal PCR swab for Chlamydia and Gonorrhoea; Trichomonas




3 Aptima PCR Urine kit for Chlamydia and Gonorrhoea



4 Universal transport media (UTM) for Herpes Simplex/Herpes Zoster, Influenza, Bordetella Pertussis and Measles PCR
(Not to be used for Chlamydia Testing)



5 Mycology pack
(for example nail clippings, skin scrapings, hair etc)



6 Surepath vial with broom and brush for cervical cytology
(High risk HPV performed from vial as appropriate)

LABELLING

Ensure all specimens are labelled with the patient's first name, last name, date of birth and date/time of collection